

# SANIGONE

SANITISE AND ITS GONE

www.sanigone.co.uk info@sanigone.co.uk















WHY SANIGONE?

## SANITISATION YOU CAN TRUST

Health professionals tell us that sanitising is the key to staying safe from Coronavirus. That is exactly what **SANIGONE** does.

Sanigone is a highly effective sanitiser and disinfectant that eliminates 99.99% of **ALL** enveloped coronaviruses together with other viruses, bacteria, yeast, funghi and mould for up to 7 days after each application.

Both of our products, **THE TOTAL RELEASE** and **THE TOTAL BLAST** also leave your space smelling clean and fresh.

**SANIGONE** is certified with three **British Standards** of testing and authentication:

BS EN 14476, BS EN 1276 and BS EN 1650.

Full certification is provided **HERE**.



**SANIGONE** is proud to be supporting a UK manufacturer especially in these difficult economic times.



**SANIGONE** cans are recyclable including the lids so cutting down on plastic waste.





## TOTAL RELEASE SANITISER

Our revolutionary **TOTAL RELEASE** room disinfectant sanitiser is a 200ml aerosol which works by dispersing it's entire contents in one application.

It is reassuringly quick and easy to use and gives immediate protection and peace of mind to everyone that the space is clean and safe.

Each aerosol covers an area of 100m<sup>3</sup> and lasts for up to 7 days giving you complete peace of mind. Even hard to reach places are sanitised during the process of application ensuring that you are protected.

#### HOW TO USE

- 1. The **SANIGONE** aerosol should be shaken well before use.
- 2. Close all doors and windows.
- 3. Ensure that the area is free of people, animals and food.
- 4. Place the aerosol on the floor in the middle of the room to be sanitised. It is best to rest it on a small dry cloth in case of any overspill.
- 5. Ensure that there are no smoke alarms directly above the aerosol.
- 6. Press and lock the total release valve located on the top of the aerosol to activate.
- 7. Exit the area and close the door.
- 8. Avoid re-entering for at least 20 minutes and ventilate if required on re entry.
- 9. The high alcohol content will ensure complete evaporation leaving no sticky residue.
- 10. NB: Once the aerosol has been activated it cannot be stopped

VIEW PRODUCT

VIEW VIDEO





## TOTAL BLAST SANITISER

Our revolutionary **TOTAL BLAST** disinfectant sanitiser spray is a 500ml aerosol which allows the user to disperse short powerful blasts of its content in to the space, onto surfaces, items and touch points that require sanitisation.

Each two to three second blast covers an area of 20cm<sup>3</sup> and ensure protection for up to 7 days giving you complete peace of mind.

Our Total Blast is perfect for sanitising areas that are more difficult to reach and can also be used as a top up for the Total Release to ensure maximum protection.

### HOW TO USE

- 1. Shake well before use
- 2. Ensure the directional spray is pointed away from the user and any other people or animals in the room
- 3. Press the nozzle for approximately 2-3 seconds to release a blast of the sanitiser which will cover an area of 20m<sup>3</sup>
- 4. As well as eliminating 99.99% of all enveloped coronaviruses, other viruses, bacteria, yeast, funghi and mould the spray delicately scents your space with pine.

VIEW PRODUCT

VIEW VIDEO





www.sanigone.co.uk

info@sanigone.co.uk